Antech Online Lab Results

Import Results to Veterinary Practice Management Software (VPMS)

Lab results can be imported directly into your veterinary practice management software (VPMS), if your VPMS vendor offers this feature. This can be done for individual results, or as a group of selected results.

From Zoasis, you will create a file that your VPMS will then read. The data from the file that you create gets plugged directly into your patient record. There are two different types of files that a VPMS could read – RTF or XML. Consult your VPMS vendor to determine which type of file you should create.

NOTE: The import feature is not available with all VPMS, and may work differently for each software system. Consult your VPMS vendor for information on availability or to troubleshoot imported records.

To Create an XML File for Import

From the main Results page, click on the arrow in the "Actions..." drop-down box and select "Import". All selected (\Box) items will be imported.

Results for 09/25/03 - 09/27/03 (new results in bold)			Page 1	Page 1 of 60 (1-10 of 599) View 💌			?
Reported	Received	<u>Patient</u>	<u>Owner</u>	<u>Doctor</u>	<u>Test</u>	Actions Print	ect
09/26/03	01/09/02	Charlie	Simmons Jean	Bruyette	DEM0942494	PDF	1
					<u>Superchem - Final</u> <u>Complete Blood Count - Fina</u> <u>T4 (RIA) - Final</u> Quantitative Platelet - Final	Save E-Mail Import	

From within an individual lab result, click on the arrow in the "Actions..." dropdown box and select "Import". Only this individual result will be imported.

Reference No. DEI	MO942494	View		Actions
Test	Results	Adult Reference Range	L	Normal Actions Print
Superchem				Save
AST (SGOT)	29	15-66 IU/L		E-Mail Import
ALT (SGPT)	36	12-118 IU/L		Consultation ^{HS}

Selecting "Import" will open a new browser window that contains the XML code for the lab report(s). This code contains all the data that your PMS will need to plug the laboratory information into the patient record. You must save this XML file to your computer, so that your PMS can read it and translate the data. In the browser toolbar, click File, then choose Save As..., as shown below. Consult your PMS vendor for instructions about where to place the file and what to name it.

File Edit View Favorites Tools Help	11					
New 🕨 🖉 Search 📷 Favorites 🏈 Media 🧭 🖏 🍙 🗐 🔚 🐮 🏦 🖧 🖓						
Open Ctrl+O	60					
Edit	30					
Save Ctrl+5	**					
Save As						
Page Setup ding="UTF-8" ?>						
Print Ctrl+P						
Print Preview						
Send						
Import and Export						
Properties ="9610" Ext-ID="9610" Status="Active">						
Work Offline Is Angeles Vet Group						
Close						
<street>1818 S Sepulveda</street>						
<street></street>						
<city>West Los Angeles</city>						
<stateprov>CA</stateprov>						
<postalcode>90025</postalcode>						
<country>USA</country>						
<phone>310-473-2951</phone>						
<td></td>						
<pre></pre>						
<pre>- <lablocation lablo="1" lablocationid="7"></lablocation></pre>						
<pre><locationcode></locationcode></pre> <pre></pre> <pre></pre> <pre></pre> <pre></pre> <pre></pre>						
- <address></address>						
<street>17672-A Cowan Avenue</street>						

To Create an RTF File for Import

From the main **Results** page, click on the arrow in the "Actions…" drop-down box and select "Save". All selected (☑) items will be downloaded into ONE file.

Results for 09/25/03 - 09/27/03 (new results in bold)			Page 1	of 60 (1-10	of 599) View 💌	Actions 💌	?
Reported	Received	<u>Patient</u>	<u>Owner</u>	<u>Doctor</u>	Test	Actions Print ^{ic}	ct
09/26/03	01/09/02	Charlie	Simmons Jean	Bruyette	DEM0942494	PDF	
					<u>Superchem - Final</u> Complete Blood Count - Final <u>T4 (RIA) - Final</u> Quantitative Platelet - Final	Save E-Mail 冹 Import	

From within an individual lab result, click on the arrow in the drop-down box and select "Save". Only this individual accession test number will be downloaded.

Reference No. DE	MO942494	View		Actions	•
Test	Results	Adult Reference Range	L	Normal Actions Print	
Superchem				Save	
AST (SGOT)	29	15-66 IU/L		E-Mail K	\$
ALT (SGPT)	36	12-118 IU/L		Consultation	n

The report will appear in your browser as an RTF (rich text format) document. Click "File" on the toolbar, and then choose "Save As…" from the menu. You can then save the file to the desired location on your computer. The file will be called "labreport.rtf" by default. You may elect to rename the file before you save it, so that you can more easily identify it at a later time.

] F	ile Edit View Inse	rt Form	at Tools Table Go To	Favorites Help			1
	New	•	Search 🛛 🙀 Favorites	🎯 History 🎁 🗟・	· 4 2 2 4	¥?	
_	Open C	Itrl+O					
Ľ	Save	hel 1 S		3 4			<u> </u>
	Save Ar						
İ	Page Setup						
	Print C	Itrl+P					- 64
. –			ITECH D IAGNOSTICS 17672-A	A Cowan Avenue Irvini	e CA 92614 Phone:	949-752-5612	- 64
	Send	•				T -1:	- 64
1	Import and Export		ANGELES VET GROUP			Tel:	- 64
· _	Import and Export					ANTECH# 027685	- 64
1	Properties		PEECO, 011 00020			ANTEON # 021000	- 64
•			Doctor	Client	Pet Name	Received	- 64
1	Create Shortcut		NONE	GOLDBERG TERI	CHESSY	09/23/2001 12:41 AM	- 64
	Delete			Case # 82784			
1	Rename		Breed	Sex	Pet Age	Reported	
1	Properties		Domestic Short Hair	SF	12	09/23/2001 07:44 AM	
-	Work Offline		Re	sults	Reference Rar	nge Units	
	Work Ornine Class		D COUNT				
			7.6		3.5-16.0	للرا ³ /للـــــــــــــــــــــــــــــــــــ	
-	RE	ЭC	8.5		5.92-9.93	ـلىر/106	
•	H	ЭB	13.	4	9.3-15.9	a/dL	